

## CROSS

### 10G Self loading stretcher



Spencer Cross is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. 10G certified when used with the fastener Spencer S-Max.

#### Caratteristiche specifiche

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Trendelenburg antishock position
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Adjustable heights system to facilitate patient transfer procedures
Foldable loading carriage to make handling in confined spaces easier

#### Optional accessories

ST42707C	S-Max fastener <small>(Mandatory for 10g versions)</small>
ST70000A	PVC anatomic mattress QMX/777
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thor.Belt w/integr. re-winding syst.
IF01047C	Track 4-30 telescopic IV holder
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

#### Standard equipment

ST70002A	STX 702 two pieces black reflex metal belts
----------	---

#### Technical data

Length	1980 ± 10 mm
Shortened length	1660 ± 10 mm
Width	580 ± 10 mm
Weight	45 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm ± 5%
Height loading carriage wheel	From 52 to 70 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

#### Class I Med.Device according to Dir. 93/42/CEE

CROSS T SELF LOADING STRETCHER  
10G YELLOW BOARD BLACK  
FRAME



**CR60045B**

**CND Classification** V0899

**Registration number** 151557

Rev.0 (22/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)