

CROSSOVER

10g Self Loading stretcher



Spencer Crossover is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover is characterized by higher loading capacity than a standard Cross.

Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Antishock positions Trendelenburg and Fowler
Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting also the front wheels.
Brakes on rear wheels
Foldable side bars Safe bar
Adjustable heights system
Foldable loading carriage

Accessories

ST42707C	S-MAX 10g fastener for stretchers <small>(Mandatory for 10g versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

PVC anatomic mattress QMX/777
STX 702 two pieces black reflex metal belts

Technical data

Length	1980 mm
Shortened length	1660 mm
Width	580 mm
Weight	45 kg
Loading capacity	250 kg
Wheels	ø200 mm
Height loading carriage wheel	From 52 to 70 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CROSS OVER TF CHROME FRAME
BLACK BOARD WITH QMX 777 10G
CR70005C



CND Classification V0899

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it