

CROSSOVER

Self loading variable heights stretcher
with black patient board and chrome frame



The importance of a positive and calm physiological attitude of the patients has been recognized and the key to success for rescue and transport. Consequently, the devices he is involved with have to shed some of their mechanical properties to become a patient friendly device; a partner rather than just an instrument to be used in the moment of need.

Here in Spencer designing and manufacturing rescue and transport devices is our specialty, but when designing the new Crossover we gave special attention also to the relationship between patient and stretcher, patient and rescuer and last but not least between rescuer and stretcher.

Our goal was the manufacture a product with ergonomic characteristics even when presented with the most complicated physical or environmental scenario. The secret behind Crossover is the Space Frame. This frame made from stainless steel and aluminum tubes specially welded has contributed to making it a truly innovative structure: Innovative not only in its shape but also in its functions where the intrinsic technological aspects make use of space normally not taken into consideration. Great attention has been paid also to the aesthetics and the harmony of the outline which is discrete without passing unobserved.

A special aspect of the Crossover is the larger diameter of both front and rear legs governed by the Twist system that now supports the functioning also of the front wheels in the same it does with the rear ones.

Equipped with the new side-bar called Safe Bar, the Crossover stretcher is supplied with two STX702 belts and mattress QMX777.

TECHNICAL DATA

Length	1970 mm
Length when reduced	1650 ± 20 mm
Width	570 mm
Weight	45 kg
Load capacity	250 kg
Front Wheels	Ø 200 with Twist
Rear Wheels	Ø 200 mm pivoting with brake
Loading height	da 480 a 750 mm
Reclining side bars	Yes
Side bar length	680 mm
Side bar height	200 mm
Frame	Steel
Lying floor	Polyethylene
Trendelenburg	Yes
Fowler	Yes
Adjustable backrest	Yes
Certified 10G	Yes

CR70005C

CROSSOVER TF BLACK
F/CROMO W/QMX777
CERTIFIED 10G



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>