

# CROSSOVER

## 10g Self Loading stretcher



Spencer Crossover is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. The loading carriage can be folded to make movement in confined spaces easier. Crossover is characterized by higher loading capacity than a standard Cross. 10G certified when used with the fastener Spencer S-Max

### Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Antishock positions Trendelenburg
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting all 4 wheels
- Brakes on rear wheels
- Foldable side bars Safe bar
- Adjustable heights system
- Foldable loading carriage

### Standard Equipment

- PVC anatomic mattress QMX/777
- STX 702 two pieces black reflex metal belts

### Technical data

|                               |                  |
|-------------------------------|------------------|
| Length                        | 1980 ± 10 mm     |
| Shortened length              | 1660 ± 10mm      |
| Width                         | 570 ± 10mm       |
| Weight                        | 45 kg            |
| Loading capacity              | 250 kg           |
| Wheels                        | ø200 mm          |
| Height loading carriage wheel | From 47 to 55 cm |
| Frame material                | Stainless steel  |
| Side bars length              | 680 mm           |
| Side bars Height              | 200 mm           |

### Accessories

|          |   |
|----------|---|
| ST42707C | S-MAX 10g fastener for stretchers <small>(Mandatory for 10g versions)</small> |
| ST00499B | STX 499 Thorax. Belt  |
| IF01047C | Track 4-30 telescopic IV holder   |
| ST00479C | FIXO KID intended for secure fixation of various age children                 |
| ST00497B | DNA STRAP Belt w/integrated re-winding system                                 |
| ST00498B | DNA STRAP Thor.Belt w/integr. re-winding syst.                                |
| ST00491A | STX 90 Telescopic headrest for stretchers                                     |

### Class I Med.Device according to Dir. 93/42/CEE

LOW CROSSOVER T F/CHROME  
P/BLACK W/QMX777 10G CER.



### CND Classification

V0899

Rev.0 (26/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)