


ADAPTER

To advance the loading carriage



Accessory that can be installed on Spencer self loading stretchers with the aim to advance the loading carriage.

Specific features

- Safe and easy to install
- Suitable for Cross and Carrera stretchers of any loading height
- The kit does not include wheels

Technical data

Material	Steel/PE
Advance amount	75 mm

EXTENSION CORD FOR LOADING
TROLLEY
CR90011B

Rev.0 (12/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it