


TRACK 2

Right IV pole for Alto



Track 2 is an IV pole for Spencer stretcher model Alto, intended to maintain the prevalence necessary in order to guarantee the continuity of the infusion operations also during the transport with the stretcher. The device must be installed on the right side of the stretcher.

Specific characteristics

- Easy to install
- Rugged and lightweight design
- It guarantees compact size through the ability to be folded up.
- The locking rings allow a quick and effective adjustment
- It can support two bags
- Bands with locking straps to prevent the fluctuations of the bags during the transport.
- Once installed, the device is maintained below the frame of the stretcher in horizontal position by a spring system

Technical data ¹

Minimum height	350 mm
Maximum height	780 mm
Materials	Al, Nylon, PE

¹ Heights measured from the stretcher frame and with IV holders opened

Class I med. Dev. according to Dir. 93/42/CEE

TRACK 2 – RIGHT IV POLE FOR
ALTO

CR90110B

CND classification



V80

Registration identification

94511

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it