

## STAR MASK – STAR MASK KEY

Face protection for CPR and keys ring



### Star Mask

An essential protective barrier, it allows a trained rescuer to have a determined approach to CPR. A must to keep with your key ring or wallet.

### Star Mask Key in nylon

Key ring in tear proof material. An extremely small space allows you to always bring along this vital device.

## TECHNICAL DATA

### Star Mask

Weight	5 g
Material	PVC Latex Free
Valve	One Way
Use mode	Disposable

### Star Mask key

Weight	20 g
Material	Nylon

CT06070B	STAR MASK KEY 2 KEY RING NYLON BLUE W/STAR MASK
----------	---



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>