

LSK 3

Laryngoscope bag



Bag suitable for storage of laryngoscopes and related accessories.

Specific features

Strap closure

Ideal to store up to 6 blades and a handle

Technical data

Dimensions	170x60x120 ± 10mm
Weight	56 g

LSK 3 EMPTY – DISPOSABLE

LARYNGOSCOPE KIT

DG00102A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it