

IME

Intubation bag



Semirigid bag to welcome laryngoscopes, blades and any other necessary accessories for intubation procedures during emergency phases.

Specific features

- Semirigid casing with zip and carrying handle
- Refractive insert to increase visibility
- Two semirigid dividers to separate the content at best
- Elastic loops on detachable mounts for attachment of the accessories anchored to the dividers by means of strap inserts

Technical data

Dimensions	310x360x100±5mm
Weight	850 ± 10 g
Materials	EVA
Detachable mounts	8 of several dimensions

IME BAG FOR INTUBATION KIT (EMPTY)
MADE OF EVA (BLACK)
DG00111A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it