

Yankauer

Sterile suction cannula



Suction cannula Yankauer without vacuum control hole.

Specific features

Extremely light
Transparent to allow easy visualization of the aspired material
Two bending zones make the access easier
Atraumatic rounded end
Supplied in sterile package

Technical data

Overall length	260 mm
Curved parts length	40mm and 50 mm
Straight part length	55 mm
Suction hole	ø 4 mm
External diameter	8 mm
Weight	About 20g

Class IIa Med. Dev. according to dir. 93/42/CEE

YANKAUER SUCTION HANDLE
WITHOUT VENT STERILE
DG00150A



CND Classification	A06010103
Registration number	335476

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it