Yankauer Vacuum Control Rigid Cannula with vacuum control





Disposable cannula equipped with on/off switch to start or stop suction. Once connected to the suction tube, it can be used with any suction device.

Specific features

Extremely light

- Transparent to allow easy visualization of the aspired material Non-slip handle
- Supplied sterile individually packaged

Technical data	
Dimensions	290x15 mm
Length curved part	60 mm
Terminal part diameter	4,5 mm
Weight	30 g
Sterilized by	EO

Class IIa Med. Dev. according to dir. 93/42/CEE	
YANKAUER – STERILE RIGID CANNULA	
W/VACUUM CONTROL	
DG00151A	0197
CND Classification	A06010103
Registration number	335477

Distributed by Spencer Italia S.r.l.

Rev.1 (23/06/2017) UNCHECKED COPY – further revisions will be available on http://support.spencer.it Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it

© Copyright - This document is property of SPENCER ITALIA SRL and can't be reproduced, used and disclosed without the written approval of Spencer Italia S.r.l.