

DGX 50

Suction catheter



Suction catheter with control valve.

Specific feature

Control valve to stop easily the suction
Size easy to identify thanks to color coded connection part

Technical data

Size	12F (4mm)
Length	550 mm
Connection color	White

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 50 – SUCTION CATHETER CH 12
DG00212C



CND Classification	R05010102
Registration number	334368

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it