

DGX 50 Suction catheter





Suction catheter with control valve..

Specific features

- Control valve to stop easily the suction
- Size easy to identify thanks to color coded connection part

Technical data	
Size	18F (6 mm)
Length	550 mm
Connection color	Red

Class IIa Med.Device according to Dir.	93/42/CEE
DGX 50 - SUCTION CATHETER CH 18	
DG00218C	LE
	0197
CND Classification	R05010102

Distributed by Spencer Italia S.r.l.

Rev.1 (10/10/2017) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u> Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

© Copyright - This document is property of SPENCER ITALIA SRL and can't be reproduced, used and disclosed without the written approval of Spencer Italia S.r.l.