

BASIC ELECTRONIC SPHYGMOMANOMETER



Electronic sphygmomanometer for non-invasive measurement of arterial blood pressure.

Specific characteristics

- Compact and fluent design
- 90 sets of memories
- Automatic power off after 3 minutes of no operation
- Low battery indication
- WHO (blood pressure classification indication)
- Systolic/Diastolic and pulse display
- Date and time display
- Irregular heart beat indicator
- Average for the latest 3 readings

Technical data

Technology	Oscillometric
Accuracy	Pressure 3 digit mmHg
Measurement range	Pressure 0 to 299 mmHg Heart rate 40 to 180 bpm
Power supply	4 AA alkaline batteries
Cuff type	Standard 22 – 32cm

Class IM MD according to Dir. 93/42/CEE

BASIC ELECTRONIC
SPHYGMOMANOMETER



DG01212E

Rev.0 (17/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it