

## DGX 410

### Aneroid sphygmomanometer



Adult sphygmomanometer for non-invasive measurement of arterial blood pressure.

#### Specific Characteristics

Measurement system with two pipes  
 The metal casing of the manometer, confers more strength  
 Bulb with adjustable screw valve  
 Manometer with numbers each 20 mmHg and symbols each 2 mmHg  
 Manometer with integrated fixing clip  
 Cuff with indications for direction of application  
 Nylon carry bag with zip closure allows rapid and easy storage

#### Technical data

Dimensions of extended cuff	530x150 mm
Suitable circumferences	From 25,4 to 40,6 cm
Tubes length	60 cm
Manometer range	0-300 mmHg
Manometer diameter	50 mm
Manometer quadrant	44 mm
Materials	Nylon, rubber, Al, Chromed brass
Weigth with bag	350 g

Class IM MD according to Dir. 93/42/CEE

DGX 410 WITH ADULT CUFF -  
 ANEROID  
 SPHYGMOMANOMETER.



**DG01410C**

<b>CND Classification</b>	C9006
<b>Registration number</b>	90110

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)