

DGX 415

Aneroid sphygmomanometer

With phonendoscope



Sphygmomanometer for non-invasive measurement of arterial blood pressure.

Specific Characteristics

Measurement system with two pipes
 The metal casing of the manometer, confers more strength
 Bulb with adjustable screw valve
 Manometer with numbers each 20 mmHg and symbols each 2 mmHg
 Manometer with integrated fixing clip
 Cuff with indications for direction of application
 Nylon carry bag with zip closure allows rapid and easy storage
 Phonendoscope with steel head and adjustable binaural arc

Technical data

| | |
|-------------------------|----------------------------------|
| Dimensions of ext. cuff | 530x150 mm |
| Suitable circumferences | From 25,4 to 40,6 cm |
| Tubes length | 60 cm |
| Manometer range | 0-300 mmHg |
| Manometer diameter | 50 mm |
| Manometer quadrant | 44 mm |
| Materials | Nylon, rubber, Al, Chromed brass |
| Phonendo head diam. | 47 mm |

Class IM MD according to Dir. 93/42/CEE

DGX 415-ANEROID
 SPHYGMOM. COMPLETE
 WITHPHONENDOSCOPE



DG01415C

CND Classification C9006

Registration number 90120

Rev.0 (31/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it