

PROFESSIONAL

Palm Aneroid sphygmomanometer



Adult sphygmomanometer for non-invasive measurement of arterial blood pressure

Specific features

- Measurement system with one pipe
- The metal casing of the manometer, confers more strength
- Bulb with adjustable screw valve
- Manometer with numbers each 20 mmHg and symbols each 2 mmHg
- Cuff with indications for direction of application
- Nylon carry bag with zip closure allows rapid and easy storage

Technical data

Dimensions of ext. cuff	530x145 mm
Suitable circumferences	From 25,4 to 40,6 cm
Tube length	50 cm
Manometer range	0-300 mmHg
Manometer diameter	68 mm
Manometer quadrant	58 mm
Materials	Nylon, PVC, Al, Chromed brass
Cuff weight	70 g
Weight sphygmomanometer	205 g
Overall weight	320 g

Class IM MD according to Dir. 93/42/CEE

PROFESSIONAL ADULT - HAND-HELD ANEROID SPHYGMOMANOMETER



DG01418A

CND Classification	C9006
Registration number	90077

Rev.2 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it