

## DGX 434

### Aneroid sphygmomanometer



Infant sphygmomanometer for non-invasive measurement of arterial blood pressure.

#### Specific Characteristics

- Measurement system with two pipes
- The metal casing of the manometer, confers more strength
- Bulb with adjustable screw valve
- Manometer with numbers each 20 mmHg and symbols each 2 mmHg
- Manometer with integrated fixing clip
- Cuff with indications for direction of application
- Nylon carry bag with zip closure allows rapid and easy storage

#### Technical data

|                    |                                  |
|--------------------|----------------------------------|
| Manometer range    | 0-300 mmHg                       |
| Manometer diameter | 50 mm                            |
| Manometer quadrant | 44 mm                            |
| Materials          | Nylon, rubber, Al, Chromed brass |

Class IM MD according to Dir. 93/42/CEE

ANEROID  
SPHYGMOMANOMETER  
INFANT



**DG01434C**

**CND Classification** C9006

**Registration number** 501169

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)