

DGX 410

Aneroid sphygmomanometer



Adult sphygmomanometer for non-invasive measurement of arterial blood pressure.

Specific Features

- Measurement system with two pipes
- The metal casing of the manometer, confers more strength
- Bulb with adjustable screw valve
- Manometer with numbers each 20 mmHg and symbols each 2 mmHg
- Manometer with integrated fixing clip
- Cuff with indications for direction of application
- Nylon carry bag with zip closure allows rapid and easy storage

Technical data

Dimensions of extended cuff	530x150 mm
Suitable circumferences	From 25,4 to 40,6 cm
Tubes length	60 cm
Manometer range	0-300 mmHg
Manometer diameter	50 mm
Manometer quadrant	44 mm
Materials	Nylon, rubber, Al, Chromed brass
Weigth with bag	350 g

Class IM MD according to Dir. 93/42/CEE

D-DGX 410 WITH ADULT CUFF -
ANEROID
SPHYGMOMANOMETER.



DG01436C

CND Classification

C9006

Distributed by Spencer Italia s.r.l.

Rev.0 (03/09/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it