

D - DGX 415

Aneroid sphygmomanometer

With phonendoscope



Sphygmomanometer for non-invasive measurement of arterial blood pressure.

Specific Characteristics

- Measurement system with two pipes
- The metal casing of the manometer, confers more strength
- Bulb with adjustable screw valve
- Manometer with numbers each 20 mmHg and symbols each 2 mmHg
- Manometer with integrated fixing clip
- Cuff with indications for direction of application
- Nylon carry bag with zip closure allows rapid and easy storage
- Phonendoscope with steel head and adjustable binaural arc

Technical data

Dimensions of ext. cuff	530x150 mm
Suitable circumferences	From 25,4 to 40,6 cm
Tubes length	60 cm
Manometer range	0-300 mmHg
Manometer diameter	50 mm
Manometer quadrant	44 mm
Materials	Nylon, rubber, Al, Chromed brass
Phonendo head diam.	47 mm

Class IM MD according to Regulation 2017/745/UE

DDGX 415-ANEROID
SPHYGMOM. COMPLETE
WITHPHONENDOSCOPE
DG01438C



Distributed by Spencer Italia s.r.l.

Rev.0 (09/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it