

## DIAGNOSTIC

### Aneroid sphygmomanometer



Adult sphygmomanometer complete of case for non-invasive measurement of arterial blood pressure.

#### Specific features

- Measurement system with two pipes
- Manometer with aluminum case for increased resistance
- Bulb with adjustable screw valve
- Nylon case with zip closure for safe and quick storage

#### Technical data

Etended cuff dimensions	530x145 mm
Case dimensions	65x180x110 mm
Suitable for circumferences	From 25,4 to 40,6 cm
Tubes length	50 cm
Manometer scale	0-300 mmHg
Dial dimensions	45 mm
Materials	Nylon, rubber, Al, Chromed brass
Weight with case	335 g

Class IM MD according to Dir. 93/42/CEE

DIAGNOSTIC ANEROIDE  
SPHYGMOMANOMETER 2  
TUBES



**DG01460D**

<b>CND Classification</b>	C9006
<b>Registration number</b>	90086

Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)