

VISO KOMPAK

Aneroid sphygmomanometer



Sphygmomanometer, for the non –invasive measurement of arterial blood pressure equipped with wall fastener.

Specific features

- Large dial
- Spiral connection pipe to sphygmomanometer to avoid obstacles while using in confined spaces
- Dial with numbering at 20mmHg and division of 2 mmHg
- Cuff with indications direction of application
- Easy to install
- Can be rotated approximately 45° left and 45° right to allow easy reading from different positions

Technical data

Extended cuff dimension	530x145 mm
Sphygmom. Dimensions	155x155x115 mm
Suitable for circumferences	from 25,4 to 40,6 cm
Pipe length	50 cm
Spiral pipe length	from 15 to 170 cm
Manometer range	0-300 mmHg
Material	PE, Nylon, PVC, Al, steel
Weight	0,68 ± 0,1 kg

Compliant with

UNI EN ISO 81060-1

MD class Im CE marked according to Dir. 93/42/CEE

VISO KOMPAK-LARGE DIAL
ANEROID SPHYGM.FOR
RAIL/WALL ATTACHMENT



DG01500A

CND classification C9006

Registration number 89887

Rev.0 (19/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it