

## VISO KOMPAK-T

Aneroid sphygmomanometer for table



Sphygmomanometer, for the non –invasive measurement of arterial blood pressure .

### Specific features

- Large dial
- The wide base makes it very stable and the inclination allows easy reading
- Dial wit numbering at 20mmHg and division of 2 mmHg
- Cuff with indications direction of application

### Technical data

Extended cuff dimension	530x145 mm
Sphygmom. Dimensions	166x176x138 mm
Suitable for circumferences	from 25,4 to 40,6 cm
Pipe length	50 cm
Manometer range	0-300 mmHg
Material	PE, Nylon, PC
Weight	0,8 ± 0,1 kg

### Compliant with

UNI EN ISO 81060-1

MD class Im CE marked according to Dir. 93/42/CEE

VISO KOMPAK-T -LARGE FACE  
ANEROID SPHYGMOM.FOR  
TABLE



**DG01502A**

**CND classification** C9006

**Registration number** 89889

Rev.1 (25/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)