

VISO BUMPER

Aneroid sphygmomanometer
with large dial and protection



Sphygmomanometer for wall installation, intended for non-invasive measurement of arterial blood pressure.

Specific features

- The dial is protected
- The connection tube to the sphygmomanometer has a spiral form to avoid obstacles while using the device in confined spaces
- Dial with numbers each 20mmHg and division of 2mmHg
- Compartment where can be placed the cuff and small accessories
- Integrated elastic band in the compartment to prevent objects come out
- The cuff has printed direction of application

Technical data

Extended cuff dimensions	530x145 mm
Viso Bumer Dimensions	ø220x150 mm
Dial diameter	ø130 mm
Circumference	From 25,4 to 40,6 cm
Tubes length	50 cm
Spiral tube length	From 15 to 170 cm
Manometer	0-300 mmHg
Materials	PE, Nylon, rubber, Al
Overall weight	1,1 ± 0,1 kg

Compliant with

UNI EN ISO 81060-1

Class IM med. Dev. According to Dir. 93/42/CEE

VISO BUMPER BLACK -
ANEROID SPHYG. WITH LARGE
DIAL AND PROTECTION



DG01601A

CND Classification	C9006
Registration number	89842

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it