

VISO BUMPER

Aneroid sphygmomanometer
with large dial and protection and 10g mounting



Sphygmomanometer for wall installation, intended for non-invasive measurement of arterial blood pressure.

Specific features

- The dial is protected
- The connection tube to the sphygmomanometer has a spiral form to avoid obstacles while using the device in confined spaces
- Dial with numbers each 20mmHg and division of 2mmHg
- Compartment where can be placed the cuff and small accessories
- Integrated elastic band in the compartment to prevent objects come out
- The cuff has printed direction of application
- The dial is protected

Compliant with

- UNI EN 1789
- UNI EN ISO 81060-1

Technical data

Extended cuff dimensions	530x145 mm
Viso Bumer Dimensions	ø220x150 mm
Dial diameter	ø130 mm
Circumference	From 25,4 to 40,6 cm
Tubes length	50 cm
Spiral tube length	From 15 to 170 cm
Manometer	0-300 mmHg
Materials	PE,Nylon, rubber, Al
Overall weight	1,1 ± 0,1 kg

Class IM med. Dev. According to Dir. 93/42/CEE

VISO BUMPER - YELLOW ANER.
SPHYGMO. LARGE DIAL 10G



DG01602A

CND Classification C9006

Registration number 89893

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it