

## VISO BUMPER

Aneroid sphygmomanometer 10g for rail



Sphygmomanometer, for the non –invasive measurement of arterial blood pressure equipped with rail fastener.

### Specific features

- Dial with protection
- Spiral connection pipe to sphygmomanometer to avoid obstacles while using in confined spaces
- Dial with numbering at 20mmHg and division of 2 mmHg
- Compartment for housing the cuff and any small accessories
- Integrated elastic band in the compartment to prevent the leakage of the accessories
- Cuff with indications direction of application
- Integrated fastener for rail bar

### Compliant with

UNI EN ISO 81060-1

### Technical data

Extended cuff dimension	530x145 mm
Viso Bumper dimension	ø220 x230 mm
Dial diameter	ø130 mm
Suitable for circumferences	from 25,4 to 40,6 cm
Pipe length	50 cm
Spiral pipe length	from 15 to 170 cm
Manometer range	0-300 mmHg
Material	PE, Nylon, rubber, Al
Weight	1,2 ± 0,1 kg

MD class Im CE marked according to Dir. 93/42/CEE

VISO BUMPER ANEROID  
SPHYGMOM. YELLOW FOR RAIL  
10G



DG01606A

CND classification C9006

Registration number 458058

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)