

VISO BUMPER

Adult cuff



Adult cuff to be used with sphygmomanometer Viso Bumper.

Specific features

- Measurement system with two pipes
- Bulb with adjustable screw valve
- Cuff with indications for direction of application

Technical data

Suitable circumferences	From 25,4 to 40,6 cm
Tubes length	50 cm
Materials	Nylon, rubber, Al, Chromed brass
Weight	170 g

VISO BUMPER – BLACK ADULT
CUFF 2 HOSES
DG01611A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it