

DGX 111

Black stethoscope



Stethoscope with 2 heads with bell and diaphragm.

Specific features

Steel head
Adjustable binaural arc

Technical data

Diaphragm diameter	47 mm
Bell diameter	30 mm
Weight	105 g
Material	Steel , PVC

Class I Med.Device according to Dir. 93/42/CEE

DGX 111 - BLACK

STETHOSCOPE

DG02310B



CND Classification C9005

Registration number 90170

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it