

DGX 601

Black cardiological phonendoscope



Phonendoscope particularly suitable for cardiac auscultation.

Specific features

- Steel head with ergonomic design
- Membranes specific for low frequencies
- Contains two pairs of soft spare eartips

Technical data

Head diameter	47 mm
Weight	170 g

Class I Med.Device according to Dir. 93/42/CEE

DGX 601 - BLACK
CARDIOLOGICAL
PHONENDOSCOPE
DG02601B



CND Classification	C9005
Registration number	90157

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it