

DGX 602

Black cardiological stethoscope



Stethoscope particularly suitable for cardiac auscultation.

Specific features

- Steel head
- Membranes specific for low frequencies
- Contains two pairs of soft spare eartips and spare diaphragm
- Protection ring on both heads to increase patient comfort

Technical data

Bell diameter	36 mm
Diaphragm diameter	47 mm
Weight	210 g

Class I Med.Device according to Dir. 93/42/CEE

VISO BUMPER - YELLOW
ANER. SPHYGMO. LARGE
DIAL 10G



DG02602B

CND Classification	C9005
Registration number	90159

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it