

L-Handle

Adult traditional light laryngoscope handle



Adult handle to be used with traditional light blades.

Specific features

Non-slip surface treatment
Needs 2 batteries type C not included

Blades connection standard

UNI EN ISO 7376

Technical data

| | |
|----------------------------|---------------|
| Materials | Chromed brass |
| Handle diameter | 30 mm |
| Overall length | 150 mm |
| Handle length | 120 mm |
| Weight (without batteries) | 195 g |

Class I Med.Dev. according to Dir. 93/42/CEE

L-HANDLE CONVENTIONAL
LIGHT ADULT
DG03100A



| | |
|---------------------|-----------|
| CND Classification | Z12021003 |
| Registration number | 141647 |

Rev.1 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it