

TFJ

Professional laryngoscope Kit



Laryngoscope kit with semi-rigid case designed for professionals in reanimation.

Specific features

- The case offers protection to the content
- Zip closure for an easy access to the content
- Laryngoscope handle with non-slip finish
- Macintosh conventional light blades
- Needs 2 batteries type C not included

Technical data

Dimensions	310x185x80 mm
Included inside the kit	
Adult laryngoscope handle	
Macintosh blades sz. 1,2,3	
Ergonomic cover for handle	
Materials	Steel, Chromed brass, EVA
Handle diameter	30 mm (without cover)
Handle length	150 mm (non-slip part 120mm)

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

TFJ LARYNG.CONV.LIGHT W/3
MACINT. BLADES SZ. 1/2/3
DG03110A



Rev.0 (20/02/2019)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it