

# LSK 1

## Laryngoscope Kit



Laryngoscope kit with semi-rigid case designed for professionals in reanimation.

### Specific features

The case offers protection to the content and is equipped with two internal detachable separators able to store blades and handles

Zip closure for an easy access to the content

Laryngoscope handle with non-slip finish

Stainless steel Macintosh and Miller blades with conventional light

Included ergonomic grip for adult handle

### Technical data

Dimensions	220x225x80 mm
Weight of the kit	1,65 ± 0,1 kg

### Included inside the kit

Adult and pediatric laryngoscope handle	
Macintosh and Miller blades sz. 0,1,2,3,4	
Ergonomic grip for adult handle	
Materials	Steel, Al, EVA

### Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 1 COMPLETE - LARYNGOSCOPE  
KIT TEN BLADES

DG03113B

Registration number



58343

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)