

LSK 2

Paediatric laryngoscope kit



A small container, created to always have at disposition an essential kit in a exceptionally reduced space made possible by meticulous organisation.
A kit, which can also be integrated with other kits which do not have a space specifically for laryngoscopes.

It is made in expansive polyethylene and Lycra, is compact and shockproof and has a zip closure. Inside there are areas for blades and handles of laryngoscopes which are well protected by the use of elastic loops.

LSK 2 contains

- a conventional light paediatric handle
- three Macintosh blades
- three Miller blades

Also available in fibre-optic version.

Dimensions: 200 x 70 x h200 mm

[8.66 x 2.75 x h8.66 in]

Weight complete: 910 g [2.00]

TECHNICAL DATA

Kit contains:

LSK 2 - PAEDIATRIC LARYNGOSCOPE KIT SIX BLADE
 LSK 2 - EMPTY PAEDIATRIC LARYNGOSCOPE
 L-HANDLE - HANDLE FOR CONVENTIONAL LIGHT
 PAEDIATRICS
 BLADE LAMA MILLER MIS.0
 S-BLADE - LAMA MILLER MIS.1
 S-BLADE - LAMA MILLER MIS.2
 C-BLADE - LAMA MACINTOSH MIS.0
 C-BLADE - LAMA MACINTOSH MIS.1
 C-BLADE - LAMA MACINTOSH MIS.2

LSK 2 complete DG03115B



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>