

LSK 4

Professional laryngoscope Kit



The image represents one size for type of blades included in the kit

Adult laryngoscope kit with semi-rigid case designed for professionals in reanimation.

Specific features

- The case offers protection to the content
- Zip closure for an easy access to the content
- Laryngoscope handle with non-slip finish
- Supplied with ergonomic silicone handle
- Macintosh and Miller conventional light blades
- Requires 2 batteries type C 1,5V not included

Technical data

Dimensions	310x185x80 mm
Weight of the kit	1440 ± 20g
Included inside the kit	
Adult laryngoscope handle	
Macintosh and Miller blades sz. 0,1,2,3,4	
Ergonomic handle	
Materials	Steel, Chromed brass, EVA
Handle diameter	30 mm (without cover)
Handle length	160 mm (non-slip part 130mm)

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 4 COMPLETE - PROFESSIONAL
LARYNGOSCOPE KIT

DG03116B

Registration number



58364

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it