

LSK 5

Holster laryngoscope Kit



Pediatric laryngoscope kit with semi-rigid case designed for professionals in reanimation

Specific features

- The case offers protection to the content
- Rear loop for fastening to the belt with strap/button closure
- Zip closure for easy access to the content
- Laryngoscope handle with non-slip finish and conventional light
- Stainless steel Macintosh blades
- Needs 2 batteries type AA not included

Technical data

Dimensions	210x145x70 mm
Weight of the kit	510 g

Content

Pediatric laryngoscope handle	
Macintosh blades sz 0,1,2	
2 pairs of non sterile latex sz L (not shown in the picture)	
Materials	Steel, Al, EVA
Handle diameter	19 mm
Length with handle	150 mm

Blade connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 5 - PAEDIATRIC
LARYNGOSCOPE KIT WITH 3
BLADES AND PAIR OF GLOVES
DG03117B



Rev.2 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it