

L-Handle

Adult fibre optic laryngoscope handle





Adult handle to be used with fibre optic blades.

Specific features
Non-slip surface treatment
High-brightness white led light
Needs 2 batteries type C not included

Chromed brass
30 mm
160 mm
130 mm
230 g

Blades	connection standard	
UNI EN	ISO 7376	

Class I Med.Dev. according to Dir. 93/42/CEE			
L-HANDLE FIBRE-OPTIC			
ADULT BLADE	((
DG03700E			
CND Classification	Z12021003		
Registration number	141661		