

L-Handle

Adult fibre optic laryngoscope handle



Adult handle to be used with fibre optic blades.

Specific features

- Non-slip surface treatment
- High-brightness white led light
- Needs 2 batteries type C not included

Blades connection standard

UNI EN ISO 7376

Technical data

Materials	Chromed brass
Handle diameter	30 mm
Overall length	160 mm
Handle length	130 mm
Weight (without batteries)	230 g

Class I Med.Dev. according to Dir. 93/42/CEE

L-HANDLE FIBRE-OPTIC
ADULT BLADE
DG03700E



CND Classification	Z12021003
Registration number	141661

Rev.1 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it