

C-Blade

Laryngoscope blade



Reusable fibre optic Macintosh blade.

Specific features

Steam sterilizable
Connection standard ISO 7376

Technical data

| | |
|--------------------|---------------------------|
| Size | 1 |
| Blade length | 68 mm |
| Total length | 93 mm |
| Weight | 76 g |
| Material | Stainless steel |
| Light transmission | Integrated optical fibers |

Class I Med.Device according to Dir. 93/42/CEE

C-BLADE - FIBRE-OPTIC
MACINTOSH BLADE SIZE 1



DG03701E

| | |
|----------------------------|-----------|
| CND Classification | Z12021003 |
| Registration number | 147180 |

Rev.0 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it