

C-Blade

Laryngoscope blade



Reusable fibre optic Macintosh blade.

Specific features

Steam sterilizable
Connection standard ISO 7376

Technical data

Size	2
Blade length	87 mm
Total length	110 mm
Weight	82 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

C-BLADE - FIBRE-OPTIC
MACINTOSH BLADE SIZE 2
DG03702E



CND Classification	Z12021003
Registration number	147181

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it