

C-Blade

Laryngoscope blade



Reusable fibre optic Macintosh blade.

Specific features

Steam sterilizable
Connection standard ISO 7376

Technical data

Size	4
Blade length	130 mm
Total length	152 mm
Weight	112 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

C-BLADE - FIBRE-OPTIC
MACINTOSH BLADE SIZE 4
DG03704E



CND Classification	Z12021003
Registration number	147183

Rev.0 (12/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it