

L-Handle

Pediatric fibre optic laryngoscope handle



Pediatric laryngoscope handle to be used with fibre optic blades.

Specific features

Non-slip surface treatment
High-brightness white led light
Needs 2 batteries type AA not included

Blades connection standard

UNI EN ISO 7376

Technical data

Materials	Alluminum
Handle diameter	19 mm
Overall length	160 mm
Handle length	130 mm
Weight (without batteries)	160 g

Class I Med.Dev. according to Dir. 93/42/CEE

L-HANDLE FIBRE-OPTIC
PAEDIATRIC
DG03705E



CND Classification	Z12021003
Registration number	141662

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it