

## L-Handle

### Pediatric fibre optic laryngoscope handle



Pediatric laryngoscope handle to be used with fibre optic blades.

#### Specific features

Non-slip surface treatment  
High-brightness white led light  
Needs 2 batteries type AA not included

#### Blades connection standard

UNI EN ISO 7376

#### Technical data

Materials	Alluminum
Handle diameter	19 mm
Overall length	160 mm
Handle length	130 mm
Weight (without batteries)	160 g

Class I Med.Dev. according to Dir. 93/42/CEE

L-HANDLE FIBRE-OPTIC  
PAEDIATRIC  
DG03705E



CND Classification	Z12021003
Registration number	141662

Rev.1 (05/12/2014)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)