

S-Blade

Laryngoscope blade



Reusable fibre optic Miller blade.

Specific features

Steam sterilizable
Connection standard ISO 7376

Technical data

Size	1
Blade length	78 mm
Total length	103 mm
Weight	66 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE FIBRE-OPTIC SIZE 1-
MILLER BLADE

DG03711E



CND Classification	Z12021003
Registration number	147198

Rev.1 (01/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it