

S-Blade

Laryngoscope blade



Reusable fibre optic Miller blade.

Specific features

Steam sterilizable
Connection standard ISO 7376

Technical data

Size	3
Blade length	172 mm
Total length	195 mm
Weight	94 g
Material	Stainless steel
Light transmission	Integrated optical fibers

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE FIBRE-OPTIC SIZE 3-
MILLER BLADE

DG03713E



CND Classification	Z12021003
Registration number	147201

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it