

D-S-Blade

Laryngoscope blade



Reusable Miller blade.

Specific features

Connection standard ISO 7376

Technical data

Size	3
Blade length	175 mm
Total length	200 mm
Weight	100 g
Material	Stainless steel
Light transmission	Traditional light

Class I MD compliant with UE Reg. 2017/745

D-S-BLADE SIZE 3-MILLER
BLADE

DG03737A

CND Classification



Z12021003

Distributed by Spencer Italia s.r.l.

Rev.0 (26/05/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it