

Professional laryngoscope Kit



Laryngoscope kit with semi-rigid case designed for professionals in reanimation.

Specific features

The case offers protection to the content
 Zip closure for an easy access to the content
 Laryngoscope handle with non-slip finish
 Macintosh fibre optic blades

Technical data

Dimensions 310x185x80 mm

Included inside the kit

Adult laryngoscope handle
 Macintosh blades sz. 2,3,4
 Ergonomic handle
 Materials Steel, Chromed brass, EVA

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Regulation 2017/745/UE

D-SET FIB.LARYNGO 1HAND., 3 FIB.
 BLADES SZ.2,3,4
 DG03740



Distributed by Spencer Italia s.r.l.

Rev.0 (25/08/2022)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it