

## LSK 4

### Laryngoscope Kit



Laryngoscope with semi-rigid case designed for professionals in reanimation.

#### Specific features

- The case offers protection to the content
- Two internal mesh pockets
- Zip closure for an easy access to the content
- Detachable separator for a better organization of the content
- Laryngoscope handle with non-slip finish
- Includes Macintosh and Miller blades
- Requires 2 batteries type C 1,5V not included

#### Technical data

Case dimensions	310x190x80 ± 10 mm
Weight of the kit	1,3 ± 0,1 Kg

#### Equipment

Fibre optic Macintosh blades sizes	1,2,3,4
Fibre optic Miller blades sizes	1,2,3,4
Materials	Steel, Chromed brass, EVA, Nylon
Handle diameter	30 mm
Handle length	160 mm (non-slip part 130mm)

#### Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 4 - PROF. LARYNGO. KIT 8  
BLADES FIBRE-OPTIC  
DG03814B



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)