

LSK 5

Laryngoscope Kit



Laryngoscope with semi-rigid case designed for professionals in reanimation.

Specific features

- The case offers protection to the content
- Rear loop for attachment to the belt with button and strap closure
- Zip closure for an easy access to the content
- Laryngoscope handle with non-slip finish
- Includes Macintosh blades
- Requires 2 batteries type C 1,5V not included

Technical data

Dimensions	210x145x70 mm
Weight of the kit	645 ± 20g

Equipment

Adult laryngoscope handle	
Macintosh blades size 1,2,3	
2 non-sterile latex gloves size L	<small>Not shown</small>
Materials	Steel, Chromed brass, EVA
Handle diameter	30 mm
Handle length	160 mm (utile 130mm)

Blades connection standard

UNI EN ISO 7376

Med.dev. CE marked according to Dir. 93/42/CEE

LSK 5 - AD.LARYNGO.KIT BLADES
FIB.OPT.W/LTX GLOVES
DG03815B



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it