

# S-Blade

## Laryngoscope blade



Reusable Miller blade.

### Specific features

Connection standard ISO 7376

### Technical data

Size	2
Material	Stainless steel
Light transmission	Traditional light

Class I Med.Device according to Dir. 93/42/CEE

S-BLADE SIZE 2-MILLER BLADE

DG04002A



CND Classification Z12021003

Registration number 147192

Rev.0 (07/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)